

PXM-20 BT Plus / PXM-40 BT Plus

The world's only hybrid powered portable X-ray

- Hybrid: Device can be operated by internal battery or by external power supply
- Great advantage for outdoor radiography
- Innovative approach by high performance lithium-ion polymer battery
- Max. 400 exposures by just one charging
- Clean diagnostic images by high frequency technology
- Strong body against external shock and contamination
- Exclusive remote controller using by hand switch



Model	PXM-20BT Plus	PXM-40BT Plus
Output		
Tube Vol./Current	100 kV / 20 mA	100 kV / 35 mA
Voltage range	40 - 100 kV in 1 kV steps	40 - 100 kV in 1 kV steps
mAs Range	0.4 mAs - 100 mAs 32 steps	0.4 mAs - 100 mAs 35 steps
Max. output	1.6 kW at 80 kV	2.4 kW at 80 kV
X-ray Tube		
Application	Toshiba D - 124	Toshiba D - 124
Focal Spot	1.2 mm	1.2 mm
Collimator		
Max. Light field: 35 cm x 35 cm @ 65 cm SID		
Lamp	LED	LED
Auto Timer	30 sec.	30 sec.
Features		
7 - segment LED read out reversible, Hand switch remote control function		
PROM Memory	15 APR	21 APR
Weight	11,2 kg	14,2 kg
Dimensions	225 x 170 x 340 mm	250 x 214 x 349 mm

PXMS-1800 Mobile Stand



- Mobile stand for portable X-ray units
- Height-adjustable arm from the bottom to max. height of 180 cm
- Easy to dis/assemble for transportation by car